

# Simulator Intravenous Arm IIa

Product No. **M50B-A**

The true-to-life IV training at budget friendly price.  
Improved injection pad allows more repetitive training.

Life-like vein palpation and realistic needle tip feeling.  
Three veins for practice: Basilic, Median antebrachial and Cephalic veins.  
Injection pads are easily replaceable.  
No needle marks remain.



## Set includes:

- 1 arm model with supporter stand
- 1 injection pad for arm model
- 1 syringe
- 1 jar of coloring powder (red)
- 1 spoon
- 1 injection bottle

Specifications are subject to change.

## Replacement parts;

11267-210 replacement injection pad

## Strap-on Venipuncture Trainer

Product No. **M87**

Simple, economical training kit for daily training.  
Vein tube makes three lines, two normal and one thin, around the center of the pad.  
Realistic tissue and vein wall resistance.  
Vein tube can be filled and refilled easily using soft plastic bottle that works as a manual pump.  
Puncture pad is replaceable.

## Set includes:

- 5 puncture pads (vein tube embedded)
- 5 supporting frames with belt
- 5 injector tubes
- 5 drain tubes
- 5 injection bottles

\*Simulated blood powder is not included in the set.

## Replacement parts

- 11326-010 replacement puncture pad (1 piece)
- 11326-020 replacement puncture pad (a set of 5)

Specifications are subject to change.

## Intradermal Injection Training Simulator

Product No. **M94**

This simulator provides training of intradermal injection to the inner forearm. A clear characteristic skin wheal is formed when proper injection procedure is followed.

Life-like touch and resistance of the injection site, close to that of a real human arm, allows pulling the skin and making punctures like a real patient. Injection pads are easily exchangeable.  
No needle trace remains.

## Set includes:

- 5 puncture pads
  - 1 forearm model
  - 1 stand
  - 1 carrying case
- packing size:  
31 x 34 x 14 H cm 3.3kg

## Replacement parts

11346-010 replacement puncture pad

Specifications are subject to change.

**KYOTO KAGAKU**

15 Kitanechoya-cho, Fushimi-ku,  
Kyoto, 612-8388, JAPAN  
Tel: +81-75-605-2510  
Fax: +81-75-605-2519  
www.kyotokagaku.com  
rw-kyoto@kyotokagaku.co.jp